

SURFLINK[®]

ACCESSORIES



Carefree connectivity

SurfLink® Mobile 2 is a hands-free cell phone transmitter, assistive listening device, media streamer and hearing aid remote all rolled into one. When paired with our 900sync™ wireless hearing aids, it's everything you'll need to talk on the phone, enjoy TV, music and so much more. You'll get clear and consistent wireless streaming anytime, anywhere.

SURFLINK® 2
MOBILE



*The one
device that
does it all.*

Enjoy hands-free phone calls



Talking on the phone has never been easier with SurfLink Mobile 2.

When synced with your cell phone, SurfLink Mobile 2 can activate incoming calls, control volume and more. In fact, your phone never needs to leave your pocket. Your hearing aids can act as both microphone (picking up your voice and sending it to the person you're talking to) and receiver (their voice streams directly to your hearing aids).





Enhance one-on-one conversations



Use SurfLink Mobile 2's powerful microphone to better isolate and pick up your companion's voice in noisy environments and transmit it directly to your wireless hearing aids. Going out to eat just got fun again!



*Zeroes in on
conversations in
noisy places.*

Makes a great conference microphone



SurfLink Mobile 2 is designed to enhance listening in group settings and meetings. Simply place it in the middle of the table and the speech from those around you goes directly to your hearing aids — making it easier to hear what everyone around the table is saying!

*Picks up each
speaker in a meeting
environment.*





A convenient assistive listening device



With SurfLink Mobile 2, you can hear presenters and teachers almost as if you were next to them. Just place SurfLink Mobile 2 on the podium, or ask the speaker to wear it using the optional lanyard or clip, and it transmits their voice directly to your hearing aids. You won't miss a thing!

*Transmits the speaker's
voice directly to your
hearing aids.*

Listen to TV at your own volume

*Streams audio directly
from a TV to your
hearing aids.*



SurfLink Mobile 2 lets you stream audio from a TV directly to your

wireless hearing aids in three ways:

1. Wireless Bluetooth® connection
2. Audio-out connector cable
3. By simply placing the directional microphone near the TV speaker.

Now, you don't have to turn the volume up loud in the room to hear it at a comfortable level.

You'll be happy — and so will your companions!





SurfLink Mobile 2 is music to your ears



Can also connect to audio devices
via optional adapter cable



Get up close and personal to your favorite music or radio program with SurfLink Mobile 2. This powerful, portable transmitter easily and wirelessly streams sound from Bluetooth-enabled music players directly to your hearing aids — for a sound that's immersive and impressive.

*Wirelessly streams
sound from devices.*

Life Without Limitations

Combined with our innovative SurfLink accessories, wireless hearing aids can also make it easier to enjoy TV, music, cell phone conversations and more — no wires attached!

SurfLink Remote

This handy remote lets you adjust and control your wireless hearing aids with the push of a button.



SurfLink Media 2

This set-and-forget media streaming solution connects to TVs, MP3 players and more, to wirelessly stream audio directly to your hearing aids. You can hear the volume at the level you want, while everyone else in the room can listen at the volume they want.



SurfLink Remote Microphone

A discreet, lightweight microphone worn by a conversation partner to aid patients in one-on-one conversations in difficult listening environments.



SurfLink accessories are compatible with wireless products only.

Carefree
connectivity
has arrived

SURFLINK®

Our innovative line of SurfLink accessories — combined with our wireless hearing aids — let you live life the way you want, no wires attached.



All trademarks or trade names are properties of their respective owners.

© 2016 Starkey Hearing Technologies. All Rights Reserved.
85229-007 4/16 BROC0304-03-EE-SG